

STED & STED-FCS

Dr. B. Christoffer Lagerholm
Facility Manager


Wolfson Imaging Centre - Oxford
Weatherall Institute of Molecular Medicine,
University of Oxford,
John Radcliffe Hospital

christoffer.lagerholm@imm.ox.ac.uk





Live Cell Microscopy

Observation of living cells: Non-Invasive





Light + Far-Field: non-invasive!




Far-Field Fluorescence Microscopy

Resolution: Goal




Far-Field Fluorescence Microscopy


Resolution: Goal



⇒ molecular scale



10 nm






10 μm

Live Cell Far-Field Microscopy

Fluorescence





Fluorescence-Labeling!





... from chemists

... from biology

(proteins)





Excitation



Camera-Based Far-Field Microscopy



Wide-Field Setup



Liver-Cells: **Nucleus** and **Cell-skeleton**

- Large area illuminated
- **Camera detection:**
image taken in one step

Far-Field Fluorescence Microscopy *Confocal Setup*



Liver-Cells: **Nucleus** and **Cell-skeleton**


- Small area illuminated
- **Point detection: scanning** required to construct image
- **Confinement along z** (pinhole)

Far-Field Microscopy

Resolution Limit: Diffraction Barrier

Far-Field Fluorescence Microscopy: Focussing of light

- away from surfaces – inside cells (3D)




**Size of focus / point-spread function
limited by diffraction of light!!!**

$$\Delta x = \frac{\lambda}{2n \sin \alpha}$$

Ernst Abbe 1873


Far-Field Microscopy

Resolution Limit: Diffraction Barrier




Far-Field Microscopy

Surpassing the Resolution Limit: Turning ON/OFF




Far-Field Microscopy

Surpassing the Resolution Limit: Turning ON/OFF




Far-Field Microscopy

Surpassing the Resolution Limit: Turning ON/OFF







Fluorescence Microscopy

STED Microscopy






Fluorescence Microscopy

STED Microscopy




Fluorescence Microscopy

STED Microscopy



Fluorescence Microscopy


STED Microscopy



STED Microscopy




Dynamical confinement of resolution

Nanoscale observation areas: CONTINUOUS TUNING of spatial resolution!



STED-Microscopy




Sub-Diffraction Imaging



20nm Crimson beads
633nm exc, 90ps, 30kW/cm²
785nm STED 200ps, 76MHz

STED-Microscopy




Sub-Diffraction Imaging



20nm Crimson beads
633nm exc, 90ps, 30kW/cm²
785nm STED 200ps, 76MHz

STED-Microscopy




Sub-Diffraction Imaging



20nm Crimson beads
633nm exc, 90ps, 30kW/cm²
785nm STED 200ps, 76MHz

STED-Microscopy




Sub-Diffraction Imaging



20nm Crimson beads
633nm exc, 90ps, 30kW/cm²
785nm STED 200ps, 76MHz


STED-Microscopy

Sub-Diffraction Imaging



20nm Crimson beads
633nm exc, 90ps, 30kW/cm²
785nm STED 200ps, 76MHz



STED-Microscopy Setup










Focal Volume Confinement

Focal Engineering – Local Zero

Phase mask:





STED PSF:




STED-Microscopy

Sub-Diffraction Imaging – 3D

3D STED nanoscopy I




Fluorescent 100nm Beads – multiple layers on cover glass



STED-Microscopy




Sub-Diffraction Imaging – 3D

3D STED nanoscopy II – iso STED



Scale bar: 250nm





Mitochondria in Vero cells:
outer membrane protein Tom20 (NK51, red)



Scale bar: 1µm

STED-Microscopy

Multi-Color Sub-Diffraction Imaging



3 laser lines:

Large Stokes shift dye – only one excitation laser
 Mitochondria in Vero cells:
 outer membrane protein Tom20 (NK51, red)
 matrix protein Hsp70 (Dy-485XL, green)

Schmidt et al NatMethods 2008

4 laser lines:

Synaptophysin (red, Atto647N) + syntaxin1 (green, Atto532)
 in neurons

Donnert et al BiophysLett 2006 / Meyer et al Small 2008



Fluorescence Nanoscopy

STED imaging on single NV centers

Ultrastable: apply very high STED intensities

→ Reach ultimate spatial resolutions!!!

$$\Delta r \approx \frac{\lambda}{2n \sin \alpha \sqrt{1 + I/I_{sat}}}$$






Commercial Leica SP8 gSTED Microscope at WIMM




- The LEICA SP8 is a conventional inverted laser-scanning confocal microscope that is additionally equipped for super-resolution STED imaging (resolution_{x,y} ~ 50 nm)
- This system is equipped with continuous wave (CW) lasers (@ 405, 458, 488, and 514 nm), a tuneable pulsed white laser (470-670 nm), a pulsed laser at 440 nm, and a high power CW laser at 592 nm for STED imaging

Golgi (TGN46) Dye: Alexa488 with Eva Wegel/Ian Dobbie





Count: 3348480 Min: 0
 Mean: 15.200 Max: 255
 StdDev: 33.249 Mode: 0 (1363503)





Count: 3348480 Min: 0
 Mean: 0.896 Max: 73
 StdDev: 2.529 Mode: 0 (2198878)

anti-TGN46 + 2° Ab-Alexa488





Muscle myocytes/Titin Z-disk/Alexa488 with Katja Gemlich

A. Confocal Raw





De-convolved (Huygens)





B. Gated STED (50% STED Power)

Raw




De-convolved (Huygens)




Centrosomes – with Alan Wainman/Jordan Raff - 130513





Confocal



gSTED50



gSTED100





Line profiles fit were fit to Sums of Two Spatial Gaussians; $\text{FWHM} = 2 (2 \ln[2])^{0.5} w$



Scale Bar = 1 μm

AchR Clusters— with Judy Cossins / David Beeson - 230113


Confocal Images



STED Images




Two-color STED Imaging Example




α -Tubulin (Abberior Star 440SX) / Nucler Pores (Atto 488)


Confocal 440/502 nm





gSTED 50% (440/502 nm)




De-convolved (Huygens)
gSTED 50% (440/502 nm)



Commercial STED Microscope – Leica SP8X (soon 3X)






	592 GATED/CW	660 GATED/CW	775 PULSED
Strength	GFP/YFP	Multicolor	Most established spectral range
Colocalization studies	+	++	+
Photostability	+	++	++
Live cell	++	+	(+)



STED Live Cell Microscopy

Problems




STED Microscopy

Cellular Imaging




protein-heavy subunit of **neurofilaments**
in the human neuroblastoma cell line
SH-SY5Y (retinoic acid-BDNF-
differentiated);
establishes cross-links to organize
and stabilize neurofilaments in axons



STED-Microscopy

Inside Living Cells





Living Cells:

Citrine, Endoplasmatic Reticulum (ER)

Live PtK2 cells


Hein, Willig, Hell PNAS 2008




STED Microscopy

Inside Living Cells - Dynamics

YFP-transgenic mouse
Hippocampal slice
CA1 neuron
(PNAS Nägerl et al 2008)
(BiophysJ 2011)




Live Mouse
YFP
(Science Berning et al 2012)



Live-Cell (inside)
Conventional dyes, GFP, ...
Two-Photon excitation

STED can also be combined with Fluorescence Correlation Spectroscopy (STED-FCS)




Bacia, K; Kim, SA; Schwille, P. (2006) Fluorescence cross-correlation spectroscopy in living cells. Nature Methods 3: 83 – 89.

Live Cell Nanoscopy

STED-FCS


STED-Microscopy: Tuning of observation area



STED-FCS
Determine transit time
for different sizes of observation areas
(different STED intensities)


Calculate
apparent diffusion coefficient:
 $D \sim \text{area} / \text{transit time}$

Dependencies: $D(\text{diameter})$
240nm \rightarrow 30/40nm
Varies for different diffusion modes




Lipid Plasma Membrane Dynamics

Nanoscale



Lipid Plasma Membrane Organization:

- Heterogeneous distribution...
- Interaction with proteins
- Interaction with cortical cytoskeleton



20nm


Small spatial
scales!!!!

Lipid Plasma Membrane Dynamics

Interactions on the Nanoscale: Nanodomains

Lipid rafts/nanodomains?

- (Transient) cholesterol/sphingolipid-enriched
- Dense molecular packing (ordered)
- Compartmentalize cellular processes




Pike, J. Lipid Res., Keystone meeting 2006

Problem:



- heterogeneous
- + highly dynamic
- small (<200 nm)

Cytoskeleton

- Membrane divided in compartments
- Proteins: fence/hindrance in diffusion path
- Hopping diffusion



Kusumi



Missing temporal/spatial resolution


- hardly any direct observation method
- highly debated

Lipid Plasma Membrane Dynamics

Fluorescence Recordings: Lipids


Phosphoglycerolipid:

Atto647N-phosphoethanolamine (PE)



Sphingolipid:


Atto647N-sphingomyelin (SM)




Live PtK2 cells:

physiological conditions

incorporation in plasma membrane




**BSA
complex**



Lipid Plasma Membrane Dynamics


Confocal Recordings




Lipid Plasma Membrane Dynamics

Confocal Recordings

Confocal: Limited spatial resolution !!!



**Relative large confocal observation area:
averages over details on nanoscale
cannot distinguish normal diffusion
from nanoscale hindered diffusion**



SM diffusion slightly prolonged but still normal




$\tau_d \approx 20 \text{ ms (PE)} / 30 \text{ ms (SM)}$

$(1/\alpha) \approx 1 \text{ (PE / SM)}$

**Slower normal diffusion
but no anomalous diffusion???**


Lipid Plasma Membrane Dynamics

STED Nanoscopy Measurement



No heterogeneities:

Fast diffusion → Limited temporal resolution!





STED Live Cell Spectroscopy

Single Lipid Dynamics - FCS

Confocal

STED – 40nm



STED (240 -> 40nm):

PE diffusion scales with area reduction

τ_d : 20 -> 0.6ms (35-fold)

and still normal

$(1/\alpha) \approx 1$

STED:

SM diffusion much longer than PE



τ_d : 30 -> 3ms (10-fold)

and anomalous


$(1/\alpha) \approx 1.5$

STED-FCS

Lipid Membrane Diffusion + Interactions: PE + SM



PE: free diffusion
(weak trapping)



SM: trapping

→ **Complex on molecular scale**

(proteins, lipid-shells, ...)

~10 ms, no movement during trapping

Cholesterol-assisted



(COase/ β -Cyclo-Dextrin/Zaragozic acid...)

Binding partner bound to cytoskeleton

(Latrunculin/Jasplakinolide/Nocodazole...)

Slight dependence on endogenous SM level



(Myriocin)



STED-FCS

Membrane Dynamic – Lipid Structure


Lipids



Specific lipid-protein interactions!

**STED-FCS:
New approach to study molecular interactions!**

Head Groups



**Binding partners?
2nd color!**

Interactions differ for different lipids!
(trapping strength, Coase+Latrunculin dependence)

But not on dye and label position!

Arp2/3 Is Critical for Lamellipodia and Response to Extracellular Matrix Cues but Is Dispensable for Chemotaxis

Congying Wu,^{1,4,7} Sreeja B. Asokan,^{1,4,7} Matthew E. Berginski,² Elizabeth M. Haynes,^{1,4} Norman E. Sharpless,^{3,4} Jack D. Griffith,^{4,5} Shawn M. Gomez,² and James E. Bear^{1,4,6,*}

¹Department of Cell and Developmental Biology

²Departments of Biomedical Engineering, Computer Science, and Pharmacology

³Department of Genetics

⁴Lineberger Comprehensive Cancer Center

⁵Department of Biochemistry



⁶Howard Hughes Medical Institute

University of North Carolina School of Medicine, Chapel Hill, NC 27599, USA


⁷These authors contributed equally to this work

*Correspondence: jbear@email.unc.edu



DOI 10.1016/j.cell.2011.12.034



Measuring Hop-Diffusion with STED-FCS




Measuring Hop-Diffusion with STED-FCS




Measuring Hop-Diffusion with STED-FCS






NRK cells



IA32 cells




Agreement of SPT with Qdots and STED-FCS




Fluorescence Microscopy

STED Microscopy





Acknowledgement



MPI, Göttingen

Stefan Hell

+ whole group

Many Slides

Christian Eggeling

Lipid Experiments

Veronika Mueller

Alf Honigmann

Christian Ringemann

Rebecca Medda

Birka Lalkens

Giuseppe Viccidomini

Haisen Ta

Andreas Schönle

Lipid labeling

Dr. V. Belov

S. Polyakova

Lipid labeling

Dr. G. Schwarzmann

Prof. K. Sandhoff


(Kekule-Inst. Bonn)


Hop diffusion


Christian Eggeling

Debora Machado-Andrade


Mathias Clausen









Absorption



Spontaneous Emission



Stimulated Emission





Fluorescence Correlation Spectroscopy


FCS

Fluorescence intensity over time

Low fluorescent concentration
 \Rightarrow diffusion of single-molecules
 = fluorescence bursts



Statistics in Time



Fluorescence Correlation Spectroscopy (FCS)

data acquisition - calculation of correlation function

data analysis – length and density of fluctuations

Fitting: anomalous sub-diffusion: $G(t_c) \sim 1/(1 + (t_c/\tau_d)^\alpha)$

\Rightarrow transit time τ_d (\sim mass, obs. area) = decay time
 $\sim d^2 / D$

\Rightarrow anomaly $1/\alpha$:

$(1/\alpha) = 1$: normal free diffusion

$(1/\alpha) > 1$: anomalous diffusion (e.g. trapping)