

# Detectors for imaging

## How PMT's work



**Chris Power**



Product & Application Sales Specialist

[Chris.power@zeiss.com](mailto:Chris.power@zeiss.com)

# Confocal Point Scanning Laser – Scanning - Confocal

Lasers are:

- Monochromatic so better for specific fluorescence
- Focussable to a single spot






# Widefield versus Confocal Laser – Scanning - Confocal



## Wide Field



camera





limited z-resolution  
thick sections  
(out-of-focus haze)

## Confocal



scanning




high z-resolution  
3D via sectioning  
(haze removed)

# Confocal Point Scanning

## Laser – Scanning - Confocal




A minute diaphragm, situated in a **conjugated focal plane**, prevents out of focus light from being detected.



# How PMT's work

## The basic principles




PMT's utilise 2 effects:

- **Photoelectric effect**, which is that electrons can be emitted from materials when they absorb energy from light. Discovered 1887 by Heinrich Heinz but used experimentally by Albert Einstein in 1905 leading to the 1921 Nobel prize for “The discovery of the law of the photoelectric effect”.
- **Secondary emission** – A high energy particle can induce the emission of secondary particles. Discovered in 1902 - Austin and Starke
- Both effects combined in 1934 by RCA – Harrison NJ.

# How PMT's work

## The basic components





Photons hit the **Photocathode** which emits electrons by the photoelectric effect

The photoelectrons are electrostatically accelerated and focused (and sometimes shuttered) by **focussing electrodes**.

The electrons impact the **dynode** and liberate a number of secondary electrons which in turn are electrostatically accelerated and focussed onto the next **dynode** in the chain. The max gain per dynode is typically about 25x

The secondary electrons from the last dynode are collected at the **anode** where they can be measured.

# How PMT's work Adjustments



Hamamatsu R6095 Gain vs Voltage

Typically only 2 adjustments can be made:

1. **Gain (Actually voltage)** - The ratio of secondary to primary electrons emitted at each dynode depends on the energy of the incident photons and is controlled by inter-electrode potentials.

Changing the gain in the middle ranges does not change sensitivity, it just changes amplification

2. **Offset** – A residual background current at the anode is usually always present but can be subtracted. Some modern LSM's calibrate this automatically

Also:

**Analogue gain** – Old fashioned multiplication of signal using analogue electronics which introduces noise. Phased out about 15 years ago.

**Digital gain** – Multiplication of signal after signal digitised, useful in multi array detectors.

# Key factors effecting performance of PMT's



## Factors effecting sensitivity and gain:

- Size of Photocathode = Large is easy to hit so good for NDD's but has a higher background
- Fill factor – If light missed the detector, a factor in older multi array detectors  
e.g. LSM 510 Meta only 80% versus 98% on the LSM 980
- Window material = Usually borosilicate glass, only worth changing for deep UV work.
- Photocathode compounds = (More detail on coming slide)
- Window arrangement = 'Head on' versus 'side on'. Side on designs are more robust and typically used in LSM
- Number and arrangement of dynodes
- Age = (More detail on coming slide)

## Contributors to noise:



- Dynode design
- Strong magnetic fields, or weak shielding – noise when mobile phones near
- Absorbed Helium
- High temperatures (More detail on coming slide)
- Cosmic rays
- Previous exposure to bright light



# Key factors effecting performance of PMT's Photocathode compounds



- Photocathode compounds = Different combinations of materials have different band gaps and sensitivities
- In 1935 the peak QE was 0.4% at 800nm (silver Oxide-caesium)
- In 1936 changing to caesium-antimony gave 12% at 400nm (1st commercial design)
- 1953 Hamamatsu founded
- Since 2003 Gallium Arsenide Phosphide (GaAsP) designs yield 45 to 56% QE




# Key factors effecting performance of PMT's

## Temperature



- At room temperature thermally generated electrons from the cathode dominate
- For all cathode types cooling below -25 gives no advantage
- For LSM work with a small photocathode area even working at room temperature will yield almost undetectable noise levels




# Key factors effecting performance of PMT's

## Age



- Its not so much the age as the mileage
- Illustrating long term stability (1 year) for SbCs dynodes as a function of mean anode current, under conditions of constant applied voltage and illumination.
- PMT sensitivity can be tested and checked.




## Day to day

- Protect the PMT from heavy oversaturation when in use – Some designs shut down to protect themselves
- Try to avoid very high gain levels, at least for prolonged periods
- Protect the PMT from light when not in use
- Keep room temperatures under 30 degrees C

## Usually not a big concern:

- Keep away from strong magnetic fields (be wary of phones)
- Keep away from cosmic rays (deep in a building basement or underground)
- Keep away from loose helium




PMT (photomultiplier, based on anode/cathode technology)

APD (avalanche photodetector, based on avalanche diode)


- PMT: PhotoMultiplier Tube
  - Most common used detector for confocal microscopy
  - Good light collection
  - Integration of photons -> low shot noise
  - Difficult to calibrate (though possible)
  - Broad spectrum coverage
  - Many different variants including GaAsP
  - Best dynamic range options
- APD: Avalanche PhotoDetector
  - Very fast responses
  - Typically used for photon counting
  - More narrow spectral coverage
  - Very high QE
  - Very limited count rate
- Hybrid detectors
  - Front end of a PMT, back end of an APD
  - The QE of a PMT – e.g. GaAsP
  - Very good time response for FLIM

# Sampling

PMT's allow arbitrary pixel sizes – To Match Nyquist sampling






The graph illustrates the scanning of a two-point object with the minimum number of sampling points needed to avoid a loss of resolution (spacing of sampling points 0.25 AU).



Oversampling

Correct Sampling

Undersampling



System can automatically recommend the correct pixel size knowing the magnification, field of view (Zoom) and wavelengths used.


Image Size:	134.7 $\mu\text{m}$ $\times$ 134.7 $\mu\text{m}$	Pixel Size:	0.09 $\mu\text{m}$
Frame Size	1584 px $\times$ 1584 px	Presets $\nabla$	
Sampling	1.0 x	Confocal	

# Sampling


PMT's allow arbitrary pixel formats




31x1 to 8192x9192 pixels  
Up to 100's of fps




Cardiomyocyte Cells Loaded with Fluo4  
Images and Samples Courtesy of  
Ben Prosser, UPENN



Maximum frequency of 6830 lines



Maximum frequency of 813,008  
(1.23 microseconds)



# Quantitative measurement of protein dynamics and interaction


## Fluorescence Correlation Spectroscopy




Objective Lens




Excitation Volume




### 1. Concentration




### 2. Diffusivity



### 3. Type of Movement







The quantum nature of light can be made visible in two ways:

a) by reducing the intensity down to the order of single photons

And


b) by shortening the observation time, despite high intensity.

The graph above illustrates case (b) – by cutting down the observation time, it is possible to resolve individual photons of the light flux in their irregular (statistical) succession.


The ability of PMT's to count this quickly allows photon timings to be measured, just consider jitter, delays, after pulsing and count rate

Counting method -

- Photon counting = Better for weak signals
- Analogue integration = Better for brighter signals
- Some detectors have both modes



- Repeated measurement
- Histogram with Counts/Time
- Fluorescent decay curve



### Requirements:



- Pulsed laser: single or multi-photon
- Single-photon detector (e.g. GaAsP PMT, HyD, APD)
- TCSPC electronics
- Software

# FLIM

Application example: single wavelength probe




## Environmental sensing (pH, ion, oxygen, lipid)



*Carsten Hille, Carsten Dosche,  
Potsdam University, Germany*

## FLIM-FRET


(intra- and intermolecular  
interaction)




*Courtesy of Sandra Orthaus, former  
member of Leibniz Institute for Age  
Research, Jena*

## Medicine

(anatomical structures,  
disease-related markers)




*Courtesy of D. Schweitzer,  
Augenlinik FSU Jena*




# LSM 980 with Quasar detection

Versatile GaAsP Spectral Array Overcomes Crosstalk with a Single Scan



## Spectral Unmixing to Overcome Crosstalk




# LSM 980 with Quasar detection


Versatile GaAsP Spectral Array Overcomes Crosstalk with a Single Scan



## Fluorophore Emission Spectrum





## Spectral Unmixing to Overcome Crosstalk



# Simultaneous Spectral Imaging of Living Plant Cells

Realtime Separation of 4 Fluorophores in Presence of Chlorophyll



Application Note: Spectral Imaging: a Powerful Tool for Confocal Multicolor Imaging in Living Plant Cells: [www.zeiss.com/microscopy](http://www.zeiss.com/microscopy)

# PMT's versus Cameras



	Cameras	PMTs	APD
Quantum Efficiency	60-95%	20-56%	80%
Noise sources	Read Noise – noise of reading the signal - fixed	NA	
	Dark Current – noise from heat - exposure time and temperature dependent	Very low	
	Photon Shot Noise – square route of signal - signal dependent	--> Same	
	Excess Noise Factor – EMCCD	NA	
	Clock Induced Charge – All but mainly observed in EMCCD	NA	
	Random Telegraph Noise - CMOS	NA	
Gain	1 to 10 <sup>3</sup> for EMCCD	10 <sup>6</sup> to 10 <sup>8</sup>	100 to 1000
Pixel formats	Fixed, typically 0.5 to 16 Megapixels	Freely configurable from 1 pixel to 64 Megapixel	
Frame rates large pixel format	Typically 10s to 100's	Single to 10's of frame per second	
Frame rates small pixel format	Up to Hundreds	Hundreds, thousands or more with line and pixel scanning	
Bandwidth	Mhz to Ghz	10s to 100s of Mhz	Low, 1Mhz
Dynamic range	Medium	High to very high depending on type	Low
Photon counting / FLIM?	Modulation only	TCSPC	
Fill factor	70-95%	95 to 100%	
Used in	Widefield, Spinning Disk, Structured Illumination, TIRF, Lightsheet etc...	Confocal, Multi Photon	

# PMT's versus Cameras

PMTs have much better SNR at the light levels in confocal than cameras, despite the QE difference



	Cameras	PMTs	APD
Quantum Efficiency	60-95%	20-56%	80%
Noise sources	Read Noise – noise of reading the signal - fixed	NA	
	Dark Current – noise from heat - exposure time and temperature dependent	Very low	
	Photon Shot Noise – square route of signal - signal dependent	...-> Same	
	Excess Noise Factor – EMCCD	NA	
	Clock Induced Charge – All but mainly observed in EMCCD	NA	
	Random Telegraph Noise - CMOS	NA	
Gain	1 to 10 <sup>3</sup> for EMCCD	10 <sup>6</sup> to 10 <sup>8</sup>	100 to 1000
Pixel formats	Fixed, typically 0.5 to 16 Megapixels	Freely configurable from 1 pixel to 64 Megapixel	
Frame rates large pixel format	Typically 10s to 100's	Single to 10's of frame per second	
Frame rates small pixel format	Up to Hundreds	Hundreds, thousands or more with line and pixel scanning	
Bandwidth	Mhz to Ghz	10s to 100s of Mhz	Low, 1Mhz
Dynamic range	Medium	High to very high depending on type	Low
Photon counting / FLIM?	Modulation only	TCSPC	
Fill factor	70-95%	95 to 100%	
Used in	Widefield, Spinning Disk, Structured Illumination, TIRF, Lightsheet etc...	Confocal, Multi Photon	

e.g. Image at 1024x1024 @ 1fps = 1,048,576 pixel per second

Camera read & dark noise =  $1.6e^-$  per pixel\* +  $0.15e^-/p/s = 1.75e^-$  per pixel

$1.75 \times 1,048,576 = 1,835,008$  false counts, roughly the same as shot noise of 3 photons per pixel

PMT read noise = Multi Alkali – Typically  $3,000e^-/p/s$  and GaAsP  $800e^-/p/s$


*The camera noise is 2,294x larger than a GaAsP*

\* = *This is true for full frame but is much higher for cameras run in the Mhz range and used as a point detector*



# LSM 980

## A scanhead build for efficiency




# Importance of Gentle Imaging Sensitive Imaging System



nature


## Using the NoiSee workflow to measure signal-to-noise ratios of confocal microscopes

Alexia Ferrand, Kai D. Schleicher, Nikolaus Ehrenfeuchter, Wolf Heusemann & Oliver Biehlmaier



# LSM 980 – Airyscan Detection


Superresolution, sensitivity and speed




# Airyscan

## Take advantage of spatial information

Airy pattern of  
a point-like emitter



Array detector  
of Airyscan



The offset of individual detectors to the optical axis provides **additional spatial information** in Airyscan (detectors of a „conventional“ LSM just integrate all light passing through its pinhole).

**Linear deconvolution** assigns all signals (and frequencies) recorded by individual detector elements to their appropriate locations.


### **Result:**


Isotropic 1.7-fold increase in resolving power!

*(Further reading: White paper on Airyscan)*

# Airyscan

Take advantage of spatial information







- 32 GaAsP detectors in hexagonal lattice
- Each detector approx. 0.2 AU in diameter
- Total detection area approx. 1.25 AU in diameter
- **Simultaneous improvement in resolution and signal!**




# Airyscan detection

## 1<sup>st</sup> ring of detectors



# Airyscan detection




## 2nd ring of detectors






# Airyscan detection

## 3rd ring of detectors



# Pixel reassignment

Spatial reassignment of the signal




# Airyscan processing

## Isotropic resolution improvement




170nm fluorescent beads




**Confocal microscope**  
Plan-Apochromat 63x/1.4 Oil,  
633nm illumination

Approx. resolution: 260nm



**Pixel reassignment**  
1.4x improved resolution

Approx. resolution: 185nm







**Airyscan processing**  
Up to 2 x improved resolution



Approx. resolution: 153nm

# Airyscan 2 in Superresolution Mode

## Maximum Signal-to-Noise with Simultaneous Superresolution



- Traditional Confocal @ 0.2 AU
- △ Traditional Confocal @ 1 AU
- Airyscan




Neuromuscular junction, bruchpilot, *Drosophila melanogaster*, Sample courtesy of J. Pielage, Basel, Switzerland

**Get better data faster**  
Your needs - our motivation




GaAsP - PMT



HeLa cells stained for DNA (blue), microtubules (yellow) and F-actin (magenta). A. Politi, J. Jakobi and P. Lenart, MPI for Biophysical Chemistry, Göttingen, Germany.


# Airyscan

Up to 4x improvement in SNR over GaAsP PMT @1AU





- Same sample – Stable and hard to bleach
- Identical imaging parameters other than these stated above
- All images scaled with best fit (0.4% top and bottom)

# SNR Comparison on Microtubules



# Other Options for Improving Resolution and SNR

## Pitfalls of Closing the Pinhole & Deconvolution






# Airyscan for Gentle Imaging

## Pitfalls of Closing the Pinhole & Deconvolution

4% Laser Power & DCV of 0.6 AU LSM

0.5% Laser Power & Airyscan





**8x increase in laser power to  
match Airyscan SNR**


# Scaling Airyscan 2 for Today's Model Systems with Multiplex Mode

# Airyscan 2 with Multiplex 8Y Mode




8x Parallelization for High Signal-to-Noise and Simultaneous Superresolution



# Recap LSM 880 with Airyscan Fast In 2016 Fast Mode Provided Usable Speed



Resonance Scanner – 80 FPS





Airyscan Fast Mode – 96 FPS

## Cardiomyocyte Cells with tubulin-EMTB to measure microtubule buckling

Images and Samples Courtesy of Ben Prosser, UPENN – “Detyrosinated microtubules buckle and bear load in contracting cardiomyocytes”, Science April 2016

# Multiplex Mode for Airyscan 2 Provides Larger Fields-of-View

## Maintain Resolution, SNR and Speed over Larger FOVs to Gain Context




### Cardiomyocyte Cells with tubulin-EMTB to measure microtubule buckling

Images and Samples Courtesy of Ben Prosser, UPENN – "Detyrosinated microtubules buckle and bear load in contracting cardiomyocytes", Science April 2016

# Multiplex mode for ZEISS LSM 980 with Airyscan 2

## Fast and Gentle Confocal Superresolution Imaging of Large Model Systems




Mode	Confocal	Airyscan SR	Multiplex SR-4Y	Multiplex SR-8Y	Multiplex CO-8Y
Parallelization	1	1	4	8	8
Max. resolution X/Y nm	—	120/120	140/140	120/160	Confoc. or better
Max. resolution Z nm	—	350	450	450	Confoc. or better
Max. FPS @ 512x512	6.1	4.7	25.0	47.5	34.4
FPS @ max FOV (mm)	0.4	0.2	1.0	2.0	9.6
	(@Zoom 0.6 / SF20)	(@Zoom 1.7 / SF7)	(@Zoom 1 / SF12)	(@Zoom 1 / SF12)	(@Zoom 1 / SF12)
Processing 1kx1k /150 slices	—	< 30s	< 30s	< 30s	< 30s
SNR vs convent. confocal	—	4-8x	4-8x + speed	4x + speed	—

# Fast Adaptive Deconvolution of Confocal Images with Hybrid Detectors


Lower Magnification Traditional Confocal with Closed PH Requires Slow Acquisition Rates



61 minute acquisition – 258 $\mu$ m x 258 $\mu$ m x 34 $\mu$ m




Automatic and  
adaptive  
deconvolution  
with decision  
masks





Mouse Brain Section PFA Fixation; DRAQ5 and cytoplasmic GFP

# Airyscan 2 with Multiplex Mode

Gentlest and Fastest Confocal Superresolution Imaging over large volumes




4 minute acquisition –  $467\mu\text{m} \times 467\mu\text{m} \times 34\mu\text{m}$



**>15x faster over  
3.3x larger volume  
with 2 channels**

8Y-SR Multiplex mode

61 minute acquisition –  $258\mu\text{m} \times 258\mu\text{m} \times 34\mu\text{m}$




Traditional Confocal with automatic and adaptive deconvolution with decision masks




# Airyscan 2 with Multiplex Mode


Gentlest and Fastest Confocal Superresolution Imaging over large volumes




4 minute acquisition –  $467\mu\text{m} \times 467\mu\text{m} \times 34\mu\text{m}$




61 minute acquisition –  $258\mu\text{m} \times 258\mu\text{m} \times 34\mu\text{m}$



 **50x increase in  $\mu\text{m}^3/\text{sec}$**   
Retaining all spatial information

# Airyscan 2 with Multiplex

## Parallelization – Gentle and Fast imaging







Courtesy of Sarita Patnaik, PhD, Univ. of Mainz



- Cells treated with cytotoxic drug
- Gentle superresolution imaging allows cells to recover after the noxious substance is removed
- Time lapse imaging for 20 hours every 10 min

# Detectors used in the LSM 980



# Detectors used in the LSM 980



Type	Where	Why	Confocal	NLO	FLIM	FCS
Diode detector	Transmitted light path	Cheap and very robust	X	x		
Traditional Multi Alkali	Transmission PMT & Quasar channel 1	Robust to high light, very high dynamic range, better blue sensitivity	X	x		
Multi alkali - Cooled	Quasar channels 3	Better sensitivity in red.	X	x		X
GaAsP	Quasar	Better general sensitivity	X	x		X
Multi Array Multi Alkali	Quasar	32 channels of simultaneous spectral detection, used in the 510, 710 and 880, now obsolete	X	x		X
Multi Array GaAsP	Quasar	Better general sensitivity, 6 or 32 channels of simultaneous spectral detection	X	x		X
Multi Array GaAsP	AiryScan	Super resolution & super sensitive imaging	X	x		
NDD Cooled Multi alkali	NDD Reflected light or Transmitted light	Better sensitivity in blue and red.		X		
Nosepiece GaAsP	NDD Reflected light	Best sensitivity in scattered light		X		
BiG (Binary GaAsP)	NDD Reflected light or Transmitted light	Better general sensitivity		X	X	
BiG (Binary GaAsP)	Confocal Direct couple port	Better general sensitivity	X	x	X	
GaAs or other Far red optimised PMT	Confocal Direct couple port	Best NIR detection	X		X	X
Hybrid detectors	Confocal Direct couple port	Fast FLIM response	x	x	X	x
Multi Channel Plate (MCP)	Confocal Direct couple port	Fastest FLIM response	x	x	X	
APD	Dedicated FCS Unit	Used to be used for sensitivity, now obsolete as GaAsP internal more cost effective and much higher count	X	x		X

# Further reading



A good basic introduction (15 pages) and the source of several diagrams used today  
<http://www.et-enterprises.com/files/file/Understanding-photomultipliers.pdf>




Complete handbook (323 pages)  
[http://www.hamamatsu.com/resources/pdf/etd/PMT\\_handbook\\_v3aE.pdf](http://www.hamamatsu.com/resources/pdf/etd/PMT_handbook_v3aE.pdf)



[www.zeiss.com/airyscan](http://www.zeiss.com/airyscan)  
Whitepapers on airyscan and Multiplex modes

<http://zeiss-campus.magnet.fsu.edu/referencelibrary/pdfs/ZeissConfocalPrinciples.pdf>





We make it visible.